



Product Information

BT620-NS 1.5 ML GRADUATED CENTRIFUGE STERILE TUBE

Description:

1.5ml Polypropylene Graduated Tubes with flat top caps.

Unlike polystyrene tube, Polypropylene tubes are chemical inert withstand centrifugation speed of 16,000g. Autoclavable. Tubes with flat top cap are easy to mark.

Features:

- Polypropylene tube
- 1.5ml
- DNase, RNase free
- With flat top cap
- Suitable for cell culture
- Autoclavable
- Storage RT



Technical Specifications

	BT620-NS
Max. volume	1.5ml
Graduations	Yes
Bottom shape	conical
Color	Natural
Max RCF	16,000 g
Plastic materials	Polypropylene
Sterile	Yes
DNase-free	Yes
RNase-free	Yes
Withstand Temperatures	-80°C~121°C
Dimensions	Outer Diameter: 10.7mm
	Length: 4.1mm
Qty./Pk	100X10