



## CERTIFICATE OF ANALYSIS

---

Product	2'-Deoxyuridine-5'-monophosphate disodium salt
Grade	High Purity
Product Code	DB0370
Formula	$C_9H_{11}N_2O_8PNa_2$
MW	352.1
CAS#	42155-08-8
Lot No	

---

---

Test Items	Specifications	Results
Appearance	White powder	
Purity (HPLC)	≥99%	
UV Content (on dry basis)	96% -104%	
Water Content (K.F.)	≤25.0%	
Ratio of Absorbance	A250/260 0.69-0.75 A280/260 0.35-0.41	
Solubility	clear colorless solution at 0.5mg/0.1ml in water	

---

Storage: -20°C