



---

## CERTIFICATE OF ANALYSIS

---

Product	Phenol, crystalline
Grade	Ultra Pure
Product Code	PB4112
Formula	C <sub>6</sub> H <sub>6</sub> O
MW	94.11
CAS#	108-95-2
Lot No	

---

Test Items	Specifications	Results
Appearance	White Crystals	
Clarity of Solution, 5% in water	Clear Colorless	
Assay	≥99.0%	
Freezing point	≥40.5°C	
Residue after evaporation	≤0.05%	
Heavy metals(as Pb)	≤10ppm	
Preservative	None	

---

Storage: 2-8°C. Keep dry. Warm to room temperature before opening.