



CERTIFICATE OF ANALYSIS

Product	SUPER BROTH AGAR
Product Code	SD7027
Appearance:	Powder
Lot No	

Test Items	Specifications
Conductivity (10% dilution pH soln., pH (7.5g/qs 100ml water)@25C	REPORT
Solubility (7.5g/qs 100ml water)	REPORT
Supports Growth of Bacteria	PASS
	PASS

Component (g/L):

Tryptone-B	35.0g
Yeast extract	20.0g
NaCl	5.0g
Agar	15.0g

Usage: 75.0g/L

Storage: Store at room temperature. Protect from moisture.